


**10 Minute DNA Release Kit-2**  
 (30 Minutes Pure DNA kit) (Cat# JZ-002)

**Description**

This Kit was designed by taking 30 minutes processing to obtain pure DNA from mice tail or any tissue. The purified DNA (260/280=1.7-2.2) is widely used for **Southern blot or PCR on the genotyping of transgenic mice**, and also used for DNA Methylation, Gene Chip, DNA Arrays, etc.

**At present, there are not any products like this Kit by taking 30 min. to obtain purified DNA from mice tail or any tissue.** The purified DNA Extract works well for transgenic mice of the single and multi-gene and to have good expression easily (**Fig. below**) in PCR reaction.

Genotyping the P53 transgenic mouse  
 With 10 Minutes DNA Release Kit-2



The expression of P53 **+/+** at number 1, 4;  
 P53 **+/-** at number 2, 3, 5, 6, 7, 8, 9, 11;  
 P53 **-/-** at number 10;



**Kit-2 work Tools:**

- Thermo-machine or water bath which can be set up the temperature to be 90-95° C;
- Centrifuge with 10,000 – 15,000g setting up at room temperature;
- 1.5ml eppendorf tubes, tips, and micro (10-200ul scalar) pipette;
- Vortex or mixture machine that is optional (referring to process in Step #3);
- Below reagents that is not included in the Kit, that are provided by yourself
  - Isopropanol (2-Propanol), molecular biology grade, minimum 99%;
  - 70 % alcohol.
  - 0.1 x TE buffer that is made by mixing 1/10 of 1x TE and 9/10 of distilled water.

**Kit-2 Content Table**

Catalog Number	DNA Release Buffer		Samples	Storage & Use	Price(\$)
	Name	Volume			
JZ-002	Kit2-B1	13.5ml	100	To kept at at 5°-29° C. Work well for 1 year from open day <b>Alarming!</b> The below 5° C is not allowed	66.90/kit
	Kit2-B2	1.5 ml			
	Kit2-B3	6 ml			

### Kit-2 Protocol Table

Step	Action	Example
Step 1	<ul style="list-style-type: none"> <li>Place a mice tail 6-8 mm or any tissue*6-20mg into an eppendorf tube with <b>Kit2-B1,135ul</b> and <b>Kit2-B2,15ul</b> , mix it well and to be capped. * it is better to cut the tail or tissue to be a few piece.</li> </ul>	<p><b>Ex. #1</b></p>  <p><b>Ex. #2</b></p> 
Step 2	<ul style="list-style-type: none"> <li>Incubate the tube at 95° C (no shaking) for 10 minutes. (15-25 minute is also allowed for older mice tail).</li> </ul>	
Step 3	<ul style="list-style-type: none"> <li>Flick the tube for 3-5 times, or vortex it 2-3 Seconds, then add <b>60ul</b> of <b>Kit2-B3</b>. Vortex it 2-3 seconds and incubate it at room temperature (RT) For 2 minutes.</li> </ul>	
Step 4	<ul style="list-style-type: none"> <li>Centrifuge it at 15,000x g for 3 min. at RT.</li> <li>Transfer supernatant <b>170-180 ul</b> to a fresh tube, then add <b>150ul</b> of <b>Isopropanol</b> (2-propanol), invert the tube for 6-9 times and to be at RT for 2 minutes,</li> </ul>	
Step 5	<ul style="list-style-type: none"> <li>Centrifuge at 15,000x g for 3 min. at RT. The DNA pellet should be seen in bottom of the tube.</li> </ul>	
Step 6	<ul style="list-style-type: none"> <li>Be careful to remove the supernatant.</li> <li>Add <b>1 ml</b> of <b>70% alcohol</b> to the tube by inverting it 4-6 times to wash the DNA pellet.</li> <li>Centrifuge the tube at 10,000x g for 2 minutes at RT.</li> </ul>	
Step 7	<ul style="list-style-type: none"> <li>Decant the 70% alcohol</li> <li>To dry the DNA pellet for 2-3 min. at RT, (see Ex #2).</li> </ul>	
Step 8	<ul style="list-style-type: none"> <li>Add <b>20-50ul</b> of <b>0.1X TE buffer</b> to dissolve the ready DNA pellet at RT For 4-5 minutes.</li> <li>Flick the tube for 6-8 times. The purified DNA would be dissolved well.</li> </ul>	
Step 9	<ul style="list-style-type: none"> <li>Measure the DNA concentration to be used for any what you want.</li> <li>It is unnecessary to measure the DNA concentration if the DNA is used for PCR only. Usually, to dilute the DNA extract 10-20 times with 0.1x TE buffer. Take 1-2ul as the templet in 20-25ul reaction volume to run PCR.</li> <li>Store the DNA at 4° C for a few weeks and at -20° C for longer using,</li> </ul>	

### Reference:

- Zhu, HJ et al: Ten Minute DNA Kits – A Novel Approach to Obtain DNA Easily In Modern Biological Science. Nature and Science 4(2): 58-70, 2006.
- Lee TY et al.; Phylogenetic analysis by RFLP and sequencing of mitochondrial DNA in a Korean Popukation. Arch Oharm Res., 2006 Jan; 29(1):88-95.
- Ren S et al. ;A Simplified Method to Prepare PCR Template DNA for Screening Of Transgenic and knockout Mice. Contemp Top Lab Anim Sci.2001: 40(2): 27-30.
- Hofstetter JR et al.; A Comparison of Recovery Methods and Tissue Sources. Biochem Mol Med. 1997; 62(2); 197-202.
- Malumbres M et al. Isolation of High Molecular Weight DNA for Reliable Genotyping of Transgenic Mice. Bio Techniques 1997; 22(6): 1114-1119.
- Anderson S et al.; Sequence and organization of human mitochondrial genome, Nature, 1981; 290(9):457-465.